


Nasal swab sample collection

- Insert the nasal swab into the nostril exhibiting the most drainage or congestion
- Gently insert the swab, until resistance is met (depending on your test, you may need to swab at greater depth)
- Rotate the swab against the nasal wall several times and hold for up to 15 seconds
- Gently remove the swab
- You may be asked to repeat the same process in the other nostril using the same swab
- Once the sample is collected, insert swab into tube and replace the cap
- Place the closed tube into the specimen bag, seal it and drop into collection box



Specific swab instructions will be provided at the test site.
